


Product Datasheet

Human Nectin-2 / CD112 Protein (orb348815)

Description	Human Nectin-2, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Gln 32 - Leu 360 (...)
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Nectin-2
Preservatives	Tris with Glycine, Arginine and NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 62.2 kDa. The protein migrates as 70-80 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_002847.1
Purity	95%
Source	Human Nectin-2, Fc Tag (orb348815) is expressed from human 293 cells (HEK293). It contains AA Gln 32 - Leu 360 (Accession # Q92692-2).
MW	62.2 kDa
Uniprot ID	Q92692
NCBI	NP_002847.1
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human Nectin-2 prot...