

Product Datasheet

CD19 antibody (orb348854)

Description	Mouse monoclonal antibody to CD19
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF
Immunogen	Human CCL (chronic lymphocytic leukemia) cells.
Target	CD19
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Application notes	This antibody was used in the identification of the B-cell epitope CD19.
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P15391
Expiration Date	12 months from date of receipt.