



Product Datasheet

OX40L antibody (orb348867)


Description	Human monoclonal antibody to OX40L
Species/Host	Human
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	Blocking, FC
Immunogen	Oxelumab binds to human OX40L, the ligand of OX40, stimulatory receptor that is expressed on activated T cells, natural killer (NK) cells and natural killer T (NKT) cells.
Target	OX40L
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Application notes	This antibody binds to the co-stimulatory OX40L, inhibiting its interaction with OX40.
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P23510
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat larg...



Immunohistochemical staining of human th...



Immunohistochemical staining of paraffin...