

Product Datasheet

SOX2 (phospho-Ser83) antibody (orb34889)

Description Rabbit polyclonal antibody to Sox2 (phospho-Ser83) is one of a set of genes that is crucial for the ...

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications DOT

Immunogen Synthetic Peptide

Target SOX2

Form/Appearance Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles

Note For research use only

Isotype Rabbit IgG

Clonality Polyclonal

MW 34310

Uniprot ID [P48431](#)

NCBI [NP_003097.1](#)

Dilution Range DOT: 1:500

Expiration Date 12 months from date of receipt.

