

Product Datasheet

CD38 antibody (orb348893)

| | |
|----------------------------|---|
| Description | Rabbit monoclonal antibody to CD38 |
| Species/Host | Mouse |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | FC, IP |
| Immunogen | Viable cells (1 x 10E7) of the human cell line, BJAB (Human Burkitt lymphoma B cell line), which were injected intraperitoneally into 8-week-old BALB/c mice. |
| Target | CD38 |
| Preservatives | PBS with 0.02% Proclin 300. |
| Concentration | 1 mg/ml |
| Storage | Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. |
| Note | For research use only |
| Application notes | A customer reports successfully using this antibody (coupled to magnetic beads) for isolation of plasma cells from the intestinal lamina propria. |
| Isotype | IgG |
| Clonality | Monoclonal |
| Purity | Purified |
| Uniprot ID | P28907 |
| Expiration Date | 12 months from date of receipt. |