

## Product Datasheet

### TNF alpha antibody (orb348900)

<b>Description</b>	Rabbit monoclonal antibody to TNF alpha
<b>Species/Host</b>	Human
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, FC
<b>Immunogen</b>	Purified recombinant human TNF (rTNF).
<b>Target</b>	TNF alpha
<b>Preservatives</b>	PBS with 0.02% Proclin 300.
<b>Concentration</b>	1 mg/ml
<b>Storage</b>	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
<b>Note</b>	For research use only
<b>Application notes</b>	This antibody can be used to detect and neutralize human TNF alpha. Unlike Adalimumab (Ab00718), Infliximab does not bind to mouse TNF alpha (Arnoult 2017, PMID: 28584008).
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Purified
<b>Uniprot ID</b>	<a href="#">P01375</a>
<b>Expiration Date</b>	12 months from date of receipt.