

Product Datasheet

CD200R antibody (orb348904)

Description	Rabbit monoclonal antibody to CD200R
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Recombinant protein - hCD200RCD4d3+4.
Target	CD200R
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	Purified
Uniprot ID	Q8TD46
Expiration Date	12 months from date of receipt.