

Product Datasheet

CD4 antibody (orb348906)

| | |
|----------------------------|--|
| Description | Rabbit monoclonal antibody to CD4 |
| Species/Host | Rat |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | FC, IHC, IHC-Fr |
| Immunogen | Rat T-cell line (NB2) transfected with the human CD4. |
| Target | CD4 |
| Preservatives | PBS with 0.02% Proclin 300. |
| Concentration | 1 mg/ml |
| Storage | Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Monoclonal |
| Purity | Purified |
| Uniprot ID | P01730 |
| Expiration Date | 12 months from date of receipt. |