

Product Datasheet

IL-2R antibody (orb348909)

Description	Rabbit monoclonal antibody to IL-2R
Species/Host	Human
Reactivity	Human, Primate
Conjugation	Unconjugated
Tested Applications	Blocking, FC
Immunogen	Cultured T-cells derived from a patient with Mycosis fungoides.
Target	IL-2R
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P01589
Expiration Date	12 months from date of receipt.