

Product Datasheet

CD4 antibody (orb348926)

Description	Rat monoclonal antibody to CD4
Species/Host	Rat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC, IHC-Fr
Immunogen	Rat T-cell line (NB2) transfected with the human CD4.
Target	CD4
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P01730
Expiration Date	12 months from date of receipt.