

## Product Datasheet

### IL4 antibody (orb348927)

<b>Description</b>	Rat monoclonal antibody to IL4
<b>Species/Host</b>	Rat
<b>Reactivity</b>	Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, FC, IHC
<b>Immunogen</b>	BSF-1 / IL-4 purified from serum-free supernatants of cultured EL-4 cells (a C57BL/6 thymoma) induced with 4 $\beta$ -phorbol 12 $\beta$ -myristate 12 $\alpha$ -acetate for 48h.
<b>Target</b>	IL-4
<b>Preservatives</b>	PBS with 0.02% Proclin 300.
<b>Concentration</b>	1 mg/ml
<b>Storage</b>	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Purified
<b>Uniprot ID</b>	<a href="#">P07750</a>
<b>Expiration Date</b>	12 months from date of receipt.