


Product Datasheet

Human Caspase 3/CASP3 ELISA Kit (orb348938)

| | |
|----------------------------|---|
| Description | Human Caspase 3/CASP3 ELISA Kit |
| Reactivity | Human |
| Range | 31.2 pg/ml - 2,000 pg/ml |
| Tested Applications | ELISA |
| Immunogen | Expression system for standard: E.coli; Immunogen sequence: S29-D175 & A183-H277 |
| Storage | Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.) |
| Note | For research use only |
| Sample Types | cell culture supernatant sand serum. |
| MW | 31608 MW |
| Uniprot ID | P42574 |
| Sensitivity | 10 pg/ml |
| Expiration Date | Please enquire. |



Human Caspase 3 PicoKine ELISA Kit stand...