


## Product Datasheet

# Mouse Neutrophil Elastase/ELA2/ELANE ELISA Kit (orb348950)

<b>Description</b>	Mouse Neutrophil Elastase/ELA2/ELANE ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	31.2 pg/ml - 2,000 pg/ml
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Expression system for standard: NS0; Immunogen sequence: S27-R260
<b>Storage</b>	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
<b>Note</b>	For research use only
<b>Application notes</b>	Dilution ratio of 1:5, concentration in serum and plasma is around 2.2 ng/ml
<b>Sample Types</b>	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
<b>MW</b>	28648 MW
<b>Uniprot ID</b>	<a href="#">Q3UP87</a>
<b>Sensitivity</b>	10 pg/ml
<b>Expiration Date</b>	Please enquire.



Mouse  
Elastase  
PicoKine  
ELISA Kit  
standa...