



## Product Datasheet

### CD71 antibody (orb348973)

<b>Description</b>	Mouse monoclonal antibody to CD71.
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	Synthetic Peptide
<b>Target</b>	TFRC
<b>Form/Appearance</b>	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Euglobin precipitation followed by dialysis against PBS.
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
<b>Note</b>	For research use only
<b>Isotype</b>	IgM
<b>Clonality</b>	Monoclonal
<b>MW</b>	84871
<b>Uniprot ID</b>	<a href="#">P02786</a>
<b>NCBI</b>	<a href="#">NP_003225.2</a> , <a href="#">NP_001121620.1</a>
<b>Dilution Range</b>	WB: 1:100-250
<b>Expiration Date</b>	12 months from date of receipt.



WB analysis of HT-1080 whole cell lysate...