

Product Datasheet

N antibody (orb349001)

Description	Rabbit polyclonal antibody to N.
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Immunogen	Synthetic Peptide
Target	N
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	46025
Uniprot ID	P59595
NCBI	NP_828858.1
Expiration Date	12 months from date of receipt.