

Product Datasheet

eNOS (phospho-Ser1177) antibody (orb34901)

Description	Rabbit polyclonal antibody to eNOS (phospho-Ser1177)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT, IF
Immunogen	Synthetic Peptide
Target	NOS3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	133275
Uniprot ID	P29474
NCBI	NP_001153581.1 , NP_001153582.1 , NP_000594.2 , NP_001153583.1
Dilution Range	DOT: 1:500, IF/ICC: 1:10-50
Expiration Date	12 months from date of receipt.