



Product Datasheet

PDX1 (phospho-Thr11) antibody (orb34905)

Description	Rabbit polyclonal to PDX1 (phospho-Thr11)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT, IF
Immunogen	Synthetic Peptide
Target	PDX1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	30771
Uniprot ID	P52945
NCBI	NP_000200.1
Dilution Range	IF/ICC: 1:100, DOT: 1:500
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of SY5Y cell...