




Product Datasheet

TRRAP antibody (orb349106)

Description	Rabbit polyclonal antibody to TRRAP.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Target	TRRAP
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	437600
Uniprot ID	Q9Y4A5
Dilution Range	IF: 1:25, WB: 1:2000, WB: 1:2000
Expiration Date	12 months from date of receipt.



WB analysis of K562 whole cell lysates (...)



Immunofluorescence analysis of Human Cer...