



Product Datasheet

CNP antibody (orb349123)

Description	Rabbit polyclonal antibody to CNP.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Target	CNP
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	47579
Uniprot ID	P09543
Dilution Range	WB: 1:2000, WB: 1:2000
Expiration Date	12 months from date of receipt.



WB analysis of human brain lysates (Lane...