




Product Datasheet

UGT2B4 antibody (orb349131)

Description	Rabbit polyclonal antibody to UGT2B4.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	UGT2B4
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	60513
Uniprot ID	P06133
Dilution Range	WB: 1:2000, WB: 1:2000, WB: 1:2000, WB: 1:2000
Expiration Date	12 months from date of receipt.



WB analysis of human liver lysate using ...



WB analysis of human cerebellum lysate u...