




Product Datasheet

DENR antibody (orb349318)

Description	Mouse monoclonal antibody to DENR.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Target	DENR
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	IgG2b, κ
Clonality	Monoclonal
MW	22092
Uniprot ID	O43583
Dilution Range	IF: 1:25, WB: 1:1000-1:2000, IHC-P: 1:25, FC: 1:25
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human sk...



Flow cytometric analysis of Hela cells
U...