




Product Datasheet

TIMP2 antibody (orb349321)

Description	Mouse monoclonal antibody to TIMP2.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Target	TIMP2
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	IgG1, κ
Clonality	Monoclonal
MW	24399
Uniprot ID	P16035
Dilution Range	IF/ICC: 1:25, FACS: 1:25, IHC-P: 1:25, WB: 1:2000
Expiration Date	12 months from date of receipt.



WB analysis of HT-1080 whole cell lysate...



Immunohistochemical staining of human ki...