

Product Datasheet

FLAG Tag antibody (orb349333)

Description	Mouse monoclonal antibody to FLAG Tag.
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Target	KLH conjugated synthetic peptide encoding HA tag (DYKDDDDKC) was used as antigen.
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	IgG1, κ
Clonality	Monoclonal
Uniprot ID	unknow
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.

