

Product Datasheet

Myc Tag antibody (HRP) (orb349403)

Description	Mouse monoclonal antibody against Myc Tag (HRP) conjugated to HRP.
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Target	This Myc Tag antibody is generated from a mouse immunized with recombinant protein.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	unknow
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.

