

## Product Datasheet

### GAPDH antibody (HRP) (orb349404)

<b>Description</b>	Mouse monoclonal antibody against GAPDH (HRP) conjugated to HRP.
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Target</b>	GAPDH
<b>Form/Appearance</b>	Purified monoclonal antibody conjugated with HRP supplied in PBS with 0.05% (W/V) Proclin 300. This antibody is purified through a protein G column, followed by size-exclusion - Membrane filtration purification.
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>MW</b>	36053
<b>Uniprot ID</b>	<b>P04406</b>
<b>Dilution Range</b>	WB: 1:40000
<b>Expiration Date</b>	12 months from date of receipt.

