



Product Datasheet

RPL14 antibody (orb349499)

Description	Rabbit polyclonal antibody to RPL14.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthesized Peptide
Target	RPL14
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	23432
Uniprot ID	P50914
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.



WB analysis of 293T/17 whole cell lysate...