


Product Datasheet

MUC1 (phospho-Thr1224) antibody (orb34959)

Description	Rabbit polyclonal to MUC1 (phospho-Thr1224)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT, IF
Immunogen	Synthetic Peptide
Target	MUC1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	122102
Uniprot ID	P15941
NCBI	NP_001018017.1 , NP_001037857.1 , NP_001037855.1 , NP_001037856.1 , NP_001191220.1 , NP_001191218.1 , NP_001191219.1 , NP_001018016.1 , NP_001037858.1 , NP_001191214.1 , NP_001191217.1 , NP_001191216.1 , NP_001191215.1
Dilution Range	DOT: 1:500, IF/ICC: 1:10-50
Expiration Date	12 months from date of receipt.



Dot blot analysis of nitrocellulose memb...