




Product Datasheet

PLIN3 antibody (orb349623)

Description	Mouse monoclonal antibody to PLIN3.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant Protein
Target	PLIN3
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	IgG2a,k
Clonality	Monoclonal
MW	47075
Uniprot ID	O60664
Dilution Range	IF: 1:25, WB: 1:2000
Expiration Date	12 months from date of receipt.



WB analysis of HeLa whole cell lysate (L...



WB analysis of HeLa whole cell lysate (L...