

Product Datasheet

NFkB p65 (phospho-Ser536) antibody (orb34981)

Description	Rabbit polyclonal antibody to NFkB p65 (phospho-Ser536)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT
Immunogen	Synthetic Peptide
Target	RELA
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	60219
Uniprot ID	Q04206
NCBI	NP_001138610.1 , NP_068810.3 , NP_001230913.1
Dilution Range	DOT: 1:500
Expiration Date	12 months from date of receipt.

