

Product Datasheet

MAPK8 (phospho-Thr183/Tyr185) antibody (orb34984)

Description	Rabbit polyclonal to MAPK8 (phospho-Thr183/Tyr185)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT
Immunogen	Synthetic Peptide
Target	MAPK8
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	48296
Uniprot ID	P45983
NCBI	NP_620637.1 , NP_620634.1 , NP_002741.1 , NP_001265476.1
Dilution Range	DOT: 1:500
Expiration Date	12 months from date of receipt.

