

Product Datasheet

PDHE1A (phospho-Ser232) antibody (orb34991)

Description	Rabbit polyclonal antibody to PDHE1A (phospho-Ser232) is the pyruvate dehydrogenase complex that cat...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT, WB
Immunogen	Synthetic Peptide
Target	PDHA1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	43296
Uniprot ID	P08559
NCBI	NP_001166926.1, NP_001166927.1, NP_000275.1, NP_001166925.1
Dilution Range	DOT: 1:500
Expiration Date	12 months from date of receipt.

