

Product Datasheet

MAPK11 (Phospho-Y182) antibody (orb350235)

Description	Rabbit polyclonal antibody to MAPK11 (Phospho-Y182)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthesized Peptide
Target	MAPK11
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	41357
Uniprot ID	Q15759
Dilution Range	WB: 1:500
Expiration Date	12 months from date of receipt.

