

Product Datasheet

SAC2 antibody (orb35030)

Description	Rabbit polyclonal antibody to SAC2. The protein activity is specific for phosphatidylinositol 4, 5-b...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	INPP5F
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	128407
Uniprot ID	Q9Y2H2
NCBI	NP_001230124.1 , NP_055752.1 , NP_001230123.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

