

Product Datasheet

CNTNAP1 antibody (FITC) (orb350375)

Description	CNTNAP1 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Contactin-associated protein 1 protein (26-356AA)
Target	CNTNAP1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P78357
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.