



## Product Datasheet

### DDX3X antibody (orb350425)


<b>Description</b>	DDX3X antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IF, IHC, IP, WB
<b>Immunogen</b>	Recombinant Human ATP-dependent RNA helicase DDX3X protein (2-662AA)
<b>Target</b>	DDX3X
<b>Preservatives</b>	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<a href="#">O00571</a>
<b>Dilution Range</b>	WB:1:500-5000, IHC-P:1:100-500, IF:1:50-500, IP:1:200-2000
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of NIH3T3(Lane 1),...



Western blot analysis of K562(Lane 1), J...



Immunofluorescence analysis of HepG2 cel...