



Product Datasheet

UBE2D1 antibody (orb350796)


Description	UBE2D1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Ubiquitin-conjugating enzyme E2 D1 protein (1-147AA)
Target	UBE2D1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P51668
Dilution Range	WB:1:200-1000, IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human brain tissue showing UBE2D1 expression.



Immunohistochemical staining of human prostate tissue showing UBE2D1 expression.



Western blot analysis of HeLa whole cell lysate showing UBE2D1 protein bands.