

Product Datasheet

CD46 antibody (orb35082)

Description Peptide affinity purified rabbit polyclonal antibody.

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications FC, IF, IHC-P, WB

Immunogen Synthetic Peptide

Target CD46

Form/Appearance Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles

Note For research use only

Isotype Rabbit IgG

Clonality Polyclonal



MW 43747

Uniprot ID [P15529](#)


NCBI [NP_758862.1](#), [NP_758863.1](#), [NP_758869.1](#), [NP_758861.1](#), [NP_002380.3](#), [NP_758871.1](#), [NP_758860.1](#), [NP_722548.1](#)

Dilution Range WB: 1:1000, IHC-P: 1:50-100, FACS: 1:10-50, IF/ICC: 1:10-50

Expiration Date 12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...



Immunohistochemical staining of paraffin...