

Product Datasheet

gag antibody (HRP) (orb350913)

Description	gag antibody (HRP)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human immunodeficiency virus type 1 group M subtype B Gag polyprotein protein (2-132AA)
Target	gag
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P04591
Expiration Date	12 months from date of receipt.