

Product Datasheet

Streptavidin antibody (HRP) (orb350917)

Description	Streptavidin antibody (HRP)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Streptomyces avidinii Streptavidin protein (25-183AA)
Target	Streptavidin
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P22629
Expiration Date	12 months from date of receipt.