

Product Datasheet

stxB antibody (FITC) (orb350926)

Description	stxB antibody (FITC)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	FITC
Immunogen	Recombinant Enterobacteria phage H19B Shiga-like toxin 1 subunit B protein (21-89AA)
Target	stxB
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P69179
Expiration Date	12 months from date of receipt.