

Product Datasheet

Major pollen allergen Jun a 1 antibody (orb350928)

Description	Major pollen allergen Jun a 1 antibody
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Juniperus ashei Pectate lyase 1 protein (22-367AA)
Target	Pectate lyase 1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P81294
Expiration Date	12 months from date of receipt.