

Product Datasheet

Major pollen allergen Jun a 1 antibody (HRP) (orb350929)

Description	Major pollen allergen Jun a 1 antibody (HRP)
Species/Host	Rabbit
Reactivity	Plant
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Juniperus ashei Pectate lyase 1 protein (22-367AA)
Target	Pectate lyase 1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P81294
Expiration Date	12 months from date of receipt.