

Product Datasheet

Allergen Api g 1.0201 antibody (orb350932)

Description	Allergen Api g 1.0201 antibody
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Apium graveolens Major allergen Api g 1, isoallergen 2 protein (1-159AA)
Target	Major allergen Api g 1, isoallergen 2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P92918
Expiration Date	12 months from date of receipt.