

Product Datasheet

gen Amb a I antibody (HRP) (orb350993)

Description	gen Amb a I antibody (HRP)
Species/Host	Rabbit
Reactivity	Plant
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Ambrosia artemisiifolia Pectate lyase 5 protein (26-396AA)
Target	Pectate lyase 5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P27759
Expiration Date	12 months from date of receipt.