

Product Datasheet

G antibody (HRP) (orb351058)

Description	G antibody (HRP)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human respiratory syncytial virus A Major surface glycoprotein G protein (67-298AA)
Target	G
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P03423
Expiration Date	12 months from date of receipt.