

## Product Datasheet

### Pectinase antibody (HRP) (orb351210)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Pectinase antibody (HRP)  |
| <b>Species/Host</b>        | Rabbit  |
| <b>Reactivity</b>          | Plant   |
| <b>Conjugation</b>         | HRP   |
| <b>Tested Applications</b> | ELISA   |
| <b>Immunogen</b>           | Recombinant Chamaecyparis obtusa Polygalacturonase protein (51-514AA) |
| <b>Target</b>              | Polygalacturonase   |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300.<br>Constituents: 50% Glycerol        |
| <b>Form/Appearance</b>     | Liquid  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.         |
| <b>Note</b>                | For research use only   |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>Uniprot ID</b>          | <a href="#">Q7M1E7</a>  |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                       |