

Product Datasheet

Tmprss4 antibody (Biotin) (orb351328)

Description	Tmprss4 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Mouse Transmembrane protease serine 4 protein (52-435AA)
Target	Tmprss4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8VCA5
Expiration Date	12 months from date of receipt.