



Product Datasheet

SLC13A4 antibody (orb351601)


Description	SLC13A4 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Solute carrier family 13 member 4 protein (134-273AA)
Target	SLC13A4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UKG4
Dilution Range	WB:1:200-1000, IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of 293T whole cell...



Immunohistochemical staining of human te...



Immunofluorescence analysis of 293 cells...