


Product Datasheet

STARD10 antibody (orb351617)

Description	STARD10 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human START domain-containing protein 10 protein (1-291AA)
Target	STARD10
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9Y365
Dilution Range	WB:1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of 293T(Lane 1), k...