

Product Datasheet

Intermediate capsid protein VP6 antibody (FITC) (orb351855)

Description	Intermediate capsid protein VP6 antibody (FITC)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	FITC
Immunogen	Recombinant Rotavirus A Intermediate capsid protein VP6 protein (1-397AA)
Target	Intermediate capsid protein VP6
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	B1NKU0
Expiration Date	12 months from date of receipt.